

**ATCC Medium: 1065 LB Agar/Broth, Miller**

**Agar Medium**

LB Agar (BD 244520).....40 g  
DI Water.....1000 ml

Autoclave at 121°C.

**Broth Medium**

LB Broth (BD 244620).....25.0 g  
DI Water.....1000 ml

Autoclave at 121°C.

\*This product is commercially available from BD.

To make medium from scratch, follow formulation below:

**LB Agar, Miller Composition**

Tryptone.....10.0 g  
Yeast Extract.....5.0 g  
Sodium Chloride.....10 g  
Agar.....15.0 g

Final pH 7.0 +/- 0.2.

\*Omit agar for broth medium.