

**ATCC Medium: 1245 YPD Medium**

**Broth Medium**

YPD Broth(BD 242820).....50 g  
DI Water.....1000 ml

Or for solid medium:

**Agar Medium**

YPD Agar(BD 242720).....65 g  
DI Water.....1000 ml

Dispense into appropriate vessel and autoclave at 121°C.

\*YPD is commercially available from BD. To make medium from scratch, follow formulation below:

Yeast Extract.....10.0 g  
Peptone.....20.0 g  
Dextrose.....20.0 g  
Agar (if needed).....20.0 g  
DI Water.....1000 ml

Adjust for final pH of 5.6. Autoclave at 121°C.